

Hand Positioning Splints

1) Padded Safety Mitt

Prevents clients from pulling out or removing catheters, I. V. tubes and other vital life lines. Also helps prevent contracture and keeps clients from scratching or injuring themselves or others. Firm polyester filling Velcro wrist strap. 100% cotton mesh-covered. Loop for tie-down conditions. Washable. White. Latex free.

6727



2) Sof-Gel™ Palm Shield

Helps prevent skin breakdown in the palm and between the fingers. Gel in finger separators and palm won't bottom out giving the product a longer life. Ideal for moderate to severe hand tightness. Helps prevent finger contractures. Large thumb piece provides excellent thumb positioning.

A812405	Left	1 Size
A812406	Right	1 Size



3) Palm Shield

- Foam splint helps increase comfort for people who have moderate to severe hand tightness by separating and positioning the fingers and thumb.
- Large thumb piece provides excellent thumb positioning.

A812401	Left	1 Size
A812402	Right	1 Size



4) Adjustable Ulnar Deviation Splint

Optional pre-cut Polyform splint insert is heated and moulded on the hand for secure wrist support. Finger tabs give ulnar-to-radial pull. **Size:** wrist circumference

Left	Right	Wrist Measurement	
A495LS	A495RS	S	Up to 14cm
A495LM	A495RM	M	14-19cm
A495LL	A495RL	L	19 and up



5) Terry Cloth and Plastic Hand Cone

- Adjustable elastic 2,5cm wide strap fits any hand size.
- Removable, washable cover.

1501	Measures 12.7 x 2.5 x 4cm
-------------	---------------------------



6) Finger Contracture Cushion

Separates fingers to help prevent complete closure of the hand and protect palms. Polyester-filled cushion; cotton terry cloth absorbs perspiration. Secures to client's hand with 6cm wide elastic strap. Washable.

1494	M	7.6 x 12.7cm	Pack of 3
149401	L	10 x 12.7cm	Pack of 3



7) Cone and Finger Separator

Firm cone-shaped anticontracture device conforms to the hand's natural shape while keeping fingers separated to prevent skin breakdown. Plush cotton terry cover absorbs perspiration. Adjustable strap fits all hand sizes. Machine washable. Blue.

1495



Hand Positioning Splints

1) Palm Protector

Holds fingers out to prevent digging into the palm of the hand. Fabricated from closed-cell foam, lined with a pile fabric. Pre-rolled and attached to simplify application.

Hand wash, air dry.

Individual pack	Pack of 3		
A812101	A8121	Left	1 size
A812201	A8122	Right	1 size

2) Palm Protector with Finger Separators

Prevents fingers from digging into palm while abducting fingers to decrease tone.

Fabricated from closed-cell foam, lined with a pile fabric. Soft finger separator strap weaves between fingers and can be cut to adjust fit.

Hand washable.

A81201	Left	1 size
A81202	Right	1 size

3) Palm Protector with Vinyl Pocket

- Has the quality of the standard palm protector but with a vinyl pocket situated on the palmar aspect of the splint. This is for the application of an Ezeform thermoplastic insert which enables the splint to be custom moulded for an individual fit.
- Fingers are kept in place by Softstrap strapping material and secured by hook and loop fastening.

Includes: Ezeform insert.

AA9297L	Left
AA9297R	Right

4) Pucci Hand Air Short (inflatable)

To prevent or stretch contractures in the hand, volar hygiene problems. This orthosis provides a light and passive stretch to a contracted hand. Can also be used for patients with volar hygiene problems.

- The inflatable bladder placed inside the cover can be adjusted individually to the patient.
- Place the orthosis on the contracted hand with the valve on the lateral side. Once secured it should be inflated.
- Before removing the orthosis the bladder needs to be deflated.

Includes: with pump and laundry bag.

55462601 (=403L)	Left	1 size
55462602 (=403R)	Right	1 size

5) Air Soft™ Resting Hand Splint

Supports the fingers, thumb and wrist in a functional and/or resting position and allows simple adjustment of the wrist and thumb. Lined with self-wicking material. The malleable aluminum frame allows for static progressive splinting.

Easily adjustable with no tools or heat. Hand washable. Latex free.

Left	Right	Size	MP Width
55462301	55462304	S	5.7 - 7cm
55462302	55462305	M	7.3 - 8.3cm
55462303	55462306	L	7.6 - 8.6cm



Rolyan®